



Manual Counting of Bacterial Cells

M. Zeder 2011

Using the ocular grid and counting frame rules.



Counting rules: Only count objects that intersect the counting frame but do not touch the rejection line.

Calculation of Total Cell Numbers

aka. Abundance or Total Cell Counts (TCC)

A_{Grid} : Area of the Grid [mm^2] (Depends on Microscope, Objective, and Optovar)

A_{Filter} : Effective Filter Area [mm^2]

N : Number of Cells per Grid (Average from multiple Grids)

V : Filtered Volume (including Sample and Fixative) [ml]

F_{Dilution} : Dilution Factor = (Sample Volume + Fixative Volume) / Sample Volume

TCC: Total Cell Count (Cells per ml) [ml^{-1}]

$$\text{TCC} = \frac{N \cdot A_{\text{Filter}} \cdot F_{\text{Dilution}}}{A_{\text{Grid}} \cdot V}$$